

SONEX

HOSPITAL AND LABORATORY EQUIPMENT DETERGENT

AUSTRALIAN MEDICAL DEVICE CLASS 1
ARTG No. 125528

BENEFITS

- ✓ In-built scale inhibitor
- ✓ Ideal for use in hard or soft water
- ✓ Phosphate free- Environmentally Friendly
- ✓ Concentrated formula
- ✓ Complies with AS 4187 and AS 4815

PERFORMANCE

- ✓ **SONEX** is specifically formulated for use in clinical and anaesthetic automatic washing machines
- ✓ **SONEX** ensures consistent excellent results, without staining or damage to delicate materials

SONEX is a highly active liquid detergent free of dye and odour. It is designed for machine cleaning of anaesthetic equipment, instruments, plastic and glass utensils in hospitals, pathology and research laboratories.

MAJOR ADVANTAGES

- Guaranteed quality cleaning:*** Use of only highly efficient ingredients ensures consistent quality cleaning, free from rejects.
- Consistent destaining:*** High buffered alkali level ensures sustained cleaning and destaining. Excellent rinsibility gives the brightest possible results with less water spotting and streaks after heavy loads.
- Better machine performance:*** In-built protection mechanism assures continuous trouble free machine operation with clean jets, less pump troubles and freedom from scale in machines.
- Excellent rinsibility:*** Sonex ensures rapid & complete rinsibility so that cleaned articles are free of detergent residues. Very little compound is needed for excellent results.

SONEX is alkaline, non-flammable and safe for use when used in accordance with the directions for use. In formulating **SONEX** particular care was taken to ensure rapid and complete rinsibility so that cleaned articles are free of detergent residues. This is particularly important whenever research is being conducted and articles such as micro titre dishes must be completely detergent-free after final rinsing. It is equally important that all medical instruments be completely free of detergent residues.



DIRECTIONS FOR USE

<i>MACHINE WASHING</i>	Lightly soiled articles	2.0 mL per litre
	Medium soiled articles	4.0 mL per litre
	Heavily contaminated articles	6.0 mL per litre

Where machines have automatic dispensing systems manufacturer's instructions should be followed.

CHEMICAL COMPOSITION AND PROPERTIES

Surfactant System:.....	Anionic
Colour:.....	Slightly hazy pale brown
Odour:.....	Odourless
pH:	>12
Dilute pH (0.5% dilution):.....	9-10
Foam Properties:.....	Controlled low foam
Solubility:.....	Completely soluble in water
Abrasives:.....	None
Flash Point:.....	None
Specific Gravity:.....	1.07
Free Alkalinity:.....	2% KOH
Amines:.....	None
Toxicity:.....	Mildly Alkali
Shelf Life:.....	2 years

SAFETY DATA *Information is listed in the Material Safety Data Sheet.*

PACK SIZE 2 x 5 Litre Wide neck Agchem, 2x5 Litre Natural, 15L

ENVIRONMENTAL STATEMENT

Whiteley Corporation manufactures products which comply with Australian and International environmental policy. Whiteley Corporation aims to ensure that every product is 100% biodegradable, whilst achieving optimal performance and benefit to customers. All Whiteley Corporation products are supplied in 100% recyclable containers. Whiteley Corporation and its staff are fully committed to ensuring that our Australian environment and those international environs in which our products are used are enhanced by the use of technology, and that these environments are kept intact for future generations.

© This document remains property of Whiteley Corporation Pty Ltd.
Alterations are not permitted without prior written authorisation from
Whiteley Corporation Pty Ltd.