

Surfex®

Low Level Disinfectant

ABOUT SURFEX®:

SURFEX® is intended for use on environmental surfaces and for non critical medical devices such as hospital beds and theatre trolleys. SURFEX® requires a 5 minute contact time at 20°C for disinfection of most bacteria, fungi or viruses. A 10-minute contact time is required for *Mycobacteria (TB)* or if there is a known outbreak of *Clostridium difficile*.

DIRECTIONS FOR USE:

Read the SURFEX® Safety Data Sheet (SDS) and Product Insert prior to use.

SURFEX® is proven to kill:

Spores:	Fungi:	Viruses:
Clostridium difficile	Candida albicans	Adenovirus
Clostridium sporogenes	Trichophyton mentagrophytes	Herpes simplex virus 1
		Influenza
		Parvovirus
		Norovirus
Bacteria:		
Acinetobacter baumannii	Mycobacterium terrae (TB)	Staphylococcus aureus (MRSA)
Carbapenem Resistant Klebsiella pneumonia (CRE/CPE)	Proteus vulgaris	Stenotrophomonas maltophilia
Clostridium difficile (vegetative)	Pseudomonas aeruginosa	Vancomycin - Resistant Enterococci (VRE)
Escherichia coli	Salmonella choleraesuis	

STORAGE: Storage: Store below 30°C out of direct sunlight.

SPILLS: For dry spills: sweep up spilt material. Contain spillage and then collect by wet-brushing and dispose of through an appropriate receptacle.

For wet spills: absorb with inert material and dispose of into a suitable waste receptacle. Wash any remnants to sewer.

FIRST AID: For advice, contact a Poisons Information Centre (Phone Australia 131 126; New Zealand 0800 764 766) or a doctor. If swallowed, do NOT induce vomiting. If in eyes, hold eyelids apart and flush the eye continuously with running water. Continue flushing until advised to stop by the Poisons Information Centre or a doctor, or for at least 15 minutes. If skin contact occurs, remove contaminated clothing and flush skin and hair with running water.

Additional information is listed in the Safety Data Sheet

DISPOSAL INFORMATION: Discard any remaining solution according to State/Territory Land Waste Management Authorities if applicable. For Emergencies and additional information please refer to the SDS.

DANGER: Causes serious eye damage. Harmful if swallowed. Wear eye protection.

SAFETY DIRECTIONS:

IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTRE or doctor/physician. Wash hands thoroughly after use.

Do not eat, drink or smoke when using this product.

IF SWALLOWED: Call a POISON CENTRE or doctor/physician if you feel unwell.

STEPS FOR USING 8.5g SACHET IN 500ML BOTTLE

BOTTLE FILL



Step 1

- Ensure correct PPE is worn.
- Refer to SDS.
- Fill clean bottle with warm water (30-40°C)



Step 2

- Open SURFEX® 8.5g Sachet.
- Pour SURFEX® powder into water.

Please note: Do not pour water into powder.



Step 3

- Replace cap and gently swirl bottle a few times to mix powder and water.



Step 4

- Refer to colour chart for product activation.
- Product can take 5-10minutes to activate depending on water temperature. Warm water is recommended.
- Please note: SURFEX® requires a 5 minute contact time at 20°C for disinfection.



Step 5


- Record time on log sheet.



Step 6

- Apply to surface for disinfection.

COLOUR CHART:



**Over time the product will change from a blue colour to clear. It is still 'ready for use', and can be used up to the 8 hour mark.
** The solution must be used within 2 hours if there is a known outbreak of Clostridium difficile or Mycobacteria (TB)*

Please note:

SURFEX® must be used within 8 hours of dilution. The solution must be used within 2 hours if there is a known outbreak of Clostridium difficile or Mycobacteria (TB)

Disposal

Discard any remaining solution according to State/Territory Land Waste Management Authorities if applicable. For Emergencies and additional information please refer to the SDS.

STEPS FOR USING 85g SACHET IN A BUCKET

BUCKET FILL



Step 1

- Ensure correct PPE is worn.
- Refer to SDS.
- Fill clean bucket with warm water (30-40°C)



Step 2

- Open SURFEX® 85g Sachet along tear-line.
- Pour SURFEX® powder into water.

Please note: Do not pour water into powder.



Step 3

- Gently swirl a few times using a clean mop to mix powder and water.



Step 4

- Refer to colour chart for product activation.
- Product can take 5-10minutes to activate depending on water temperature. Warm water is recommended.
- Please note: SURFEX® requires a 5 minute contact time at 20°C for disinfection.



Step 5


- Record time on log sheet.



Step 6

- Apply to surface for disinfection.

COLOUR CHART:



**Over time the product will change from a blue colour to clear. It is still 'ready for use', and can be used up to the 8 hour mark.
** The solution must be used within 2 hours if there is a known outbreak of Clostridium difficile or Mycobacteria (TB)*

Please note:

SURFEX® must be used within 8 hours of dilution. The solution must be used within 2 hours if there is a known outbreak of Clostridium difficile or Mycobacteria (TB)

Disposal

Discard any remaining solution according to State/Territory Land Waste Management Authorities if applicable. For Emergencies and additional information please refer to the SDS.