



**INSTRUMENT GRADE DISINFECTANT
– LOW LEVEL**

ARTG: 354053

V-Wipes Benefits:

- ✓ Class IIb Instrument Grade Disinfectant - Low Level Wipes
- ✓ pH neutral - safe to use on most surfaces and equipment
- ✓ Highly absorbent apertured fabric
- ✓ Leaves surfaces streak free
- ✓ Extremely durable
- ✓ Non-hazardous
- ✓ Does not contain phenolics or chlorine
- ✓ Non-flammable
- ✓ Friendly to use on metals (non corrosive)



Additional information is listed in the Safety Data Sheet available at www.whiteley.com.au

HOW TO CLEAN AND DISINFECT HARD SURFACES USING V-WIPES



Step 1: Personal Protective Equipment.

Ensure appropriate PPE is worn prior to commencing.



Step 2: Clean and Remove Soiling.

If the surface is visibly soiled, wipe the surface with a V-Wipe, until the surface is visibly clean. Dispose of V-Wipe after use.



Step 3: Disinfecting.

With a new V-Wipe, wipe the surface again ensuring the V-Wipe makes contact with all of the surface. Work from the outside inwards.



Step 4: Allow surface to dry.

Allow surface to dry, ensuring a minimum 10-minute drying time.

HOW TO CLEAN AND DISINFECT SOFT SURFACES (SUCH AS CARPETS) USING V-WIPES



Step 1: Personal Protective Equipment.

Ensure appropriate PPE is worn prior to commencing.



Step 2: Clean and Remove Soiling.

If the surface is visibly soiled, wipe the surface with a V-Wipe, until the surface is visibly clean. Dispose of V-Wipe after use.



Step 3: Disinfecting.

With a new V-Wipe, wipe the surface again ensuring the V-Wipe makes contact with all of the surface. Work from the outside inwards.



Step 4: Allow surface to dry.

Allow surface to dry, ensuring a minimum 10-minute drying time.



Watch the Educational DVD on 'How to clean and disinfect blood spills using Viraclean® & V-Wipes' at www.whiteley.com.au/education

