

MICROL

0.2 MICRON FILTERED 70% ETHANOL

Use as directed for rinsing and/or drying of instruments and surfaces.

BENEFITS

- ✓ Filtered to **0.2 µm during filling**
- ✓ Ethanol is suitable for rinsing
- ✓ Made under Device Quality System
- ✓ Conforms to new regulations
- ✓ Contains no dye or other materials

PERFORMANCE

- ✓ **MICROL** is the only product in this category made by a Medical Device Manufacturer
- ✓ Ethanol is recognised for knock down sanitising performance on surfaces

MICROL is now the benchmark standard of product for device rinsing, being made to conform not just with the latest professional guidelines, but also the proposed new 'Medical Device Harmonisation' requirements. **MICROL** introduces Australian purchases to this new standard for the first time where in the past this narrow product segment has been dominated by foreign owned companies.

MAJOR ADVANTAGES:

It is recommended that Medical Devices be rinsed with 70% Ethanol prior to use to assist the drying of surfaces following High Level Disinfection or reprocessing. Whiteley Industries have introduced **MICROL**, which is made under strictly controlled quality control provisions in the most modern specialty Medical Devices premises in Australia. The plant is licensed under ISO 13488 by the Therapeutic Goods Administration (no endorsement should be implied), and also under ISO 9002 by Standards Australia.

Guaranteed quality: Use of only highly effective ingredients and strict quality surveillance ensures consistent quality.

Better device performance: **MICROL** is an essential component of instrument maintenance and care.

Excellent results: Used as directed, **MICROL**'s performance is outstanding.

MICROL will provide sanitising on surfaces with a 30 second contact time at room temperature (20°C).

(Continued)

METHOD OF USE: MICROL should be used as directed for rinsing and / or drying of instruments and surfaces. Care should be taken with instruments and surfaces that can be affected by alcohol as damage can occur with prolonged or repeated contact.

ENVIRONMENTAL STATEMENT

Whiteley Industries manufactures products that comply with Australian and International environmental policy. Whiteley Industries aims to ensure that every product is 100% biodegradable, whilst achieving optimal performance and benefit to customers.

All Whiteley Industries products are supplied in 100% recyclable containers.

Whiteley Industries and its staff are fully committed to ensuring that our Australian environment and those international environs in which our products are used are enhanced by the use of technology, and that these environments are preserved for the enjoyment of future generations.

CHEMICAL COMPOSITION AND PROPERTIES

Colour (Concentrate):	Clear Liquid
Odour:	Ethanol Odour
Solubility:	Soluble
Specific Gravity:	0.880
Dangerous Goods Rating:	Class 3 Flammable Liquid
Non volatiles:	Nil
Shelf life:	2 Years

SAFETY

It is imperative that those using **MICROL** should read the Material Safety Data Sheet prior to first use.

There is no way of predicting an individual's response to cleaning products. All cleaning products have a potential to cause skin and other irritation in susceptible individuals. Direct physical contact with the skin should be avoided. If irritation does occur when using this product or any other cleaning solution then the individual should be removed from direct contact with the solution until a satisfactory safe working environment for that individual is determined.

PACK SIZE: 12 X 500mL

© This document remains property of Whiteley Corporation Pty Ltd.
Alterations are not permitted without prior written authorisation from
Whiteley Corporation Pty Ltd.

Microl - Page 2 of 2
Date of Issue – 07/03