

SPEEDY CLEAN

Hospital Grade Anionic Neutral Detergent Ready to Use Solution

AUSTRALIAN MEDICAL DEVICE CLASS 1
ARTG No. 125529

INTRODUCTION AND BACKGROUND:

SPEEDY CLEAN is a ready to use neutral detergent that gives a streak free finish on all surfaces. **SPEEDY CLEAN** uses Whiteley Medical's Emulsion Systems Technology, which enhances **SPEEDY CLEAN** ability to remove oils and dirt that are encountered in everyday cleaning.

SPEEDY CLEAN is non-toxic, 100% biodegradable and quality assured as Whiteley Medical is an accredited company under ISO 9001:2000. **SPEEDY CLEAN** is safe for use on all surfaces.

SPEEDY CLEAN is a Class 1 Medical Device and can be used to clean other Medical Device surfaces (such as Theatre trolleys etc).

DESCRIPTION:

SPEEDY CLEAN

- is a pale blue liquid with a powder fresh fragrance. It is ready to use and is conveniently packaged in a 500mL squeeze bottle.
- is safe to use.
- is non toxic.
- is safe on all surfaces.
- will not cause staining on fabrics or carpets.
- has a neutral pH.
- is friendly to metals (non corrosive)
- is fully biodegradable.
- Is a Class 1 Medical Device.

VERSATILE:

SPEEDY CLEAN has been formulated to replace most cleaning products and is totally compatible with Viraclean, hypochlorites and other hospital grade disinfectants. Do not mix with chemicals.

INFECTION CONTROL:

SPEEDY CLEAN aids in the Infection Control process because of its' single use package. Coupled with Viraclean this makes for simplification of the cleaning and disinfection process.



USES:

SPEEDY CLEAN can be used on all surfaces including the following areas:

- Nursery
- CSSD
- Intensive Care Units
- Theatre Equipment
- Theatre tables and benches
- Ward areas
- Telephones

CHEMICAL SPECIFICATIONS:

Surfactant system:.....	Blended
Colour:.....	Pale Blue
pH:.....	Neutral
Solubility:.....	100% in water
Environmental:.....	Completely biodegradable
Specific Gravity:.....	1.00
Abrasives:.....	Nil
Free Caustic:.....	Nil
Phosphates:.....	Nil
Toxicity:.....	Non Toxic
Occupational safety:.....	Safe to Use
Dangerous Goods Rating:.....	Nil
Shelf Life:.....	2 Years
Safety Equipment:.....	Gloves recommended

ENVIRONMENTAL STATEMENT:

Whiteley Corporation manufactures products that comply with Australian and International environmental policy. Whiteley Corporation aims to ensure that every product is 100% biodegradable, whilst achieving optimal performance and benefit to customers. All Whiteley Corporation products are supplied in 100% recyclable containers.

Whiteley Corporation and its staff are fully committed to ensuring that our Australian environment and those international environs in which our products are used are enhanced by the use of technology, and that these environments are preserved for the enjoyment of future generations.

© This document remains property of Whiteley Corporation Pty Ltd.
Alterations are not permitted without prior written authorisation from
Whiteley Corporation Pty Ltd.

Pack Size:
12 x 500mL