


PROXY

Automatic Endoscope Reprocessing (AER) Chemistry System*

Key Benefits:

- ✓ Australian made
- ✓ Extremely cost efficient

Insist on **Whiteley Medical** when ordering your Endoscopy products



PROXY A

Anticorrosive Additive

An additive for use with PROXY P.

ARTG No: 274095

✓ 2 x 5 Litres

Product Code: 2x5L - 160622



PROXY C+

Concentrated Pre-disinfectant Cleaner

An acidic pre-disinfectant cleaner, compatible with PROXY automatic endoscope reprocessing chemistry system.

ARTG No: 125527

✓ 2 x 5 Litres

Product Code: 2x5L - 160642



PROXY P

High Level Disinfectant & Sterilant

A 5% peracetic acid based high level disinfectant and sterilant. For use with PROXY A. ARTG No: 274095

✓ 2 x 5 Litres

Product Code: 2x5L - 160612



PROXY CLN

Multi-Enzyme Detergent

Multi-Enzyme, low foaming detergent compatible with PROXY automatic endoscope reprocessing chemistry system. ARTG No: 125529

✓ 2 x 5 Litres

Product Code: 2x5L - 160652



PROXY L

70% Ethanol Solution

70% solution of Ethanol in water.

ARTG No: 127945

✓ 4 x 2 Litres

Product Code: 4x2L - 160632

*Suitable for use in Soluscope 3 PAA AERs, please contact your Whiteley Medical representative prior to first use.

For more information about Whiteley Medical products visit www.whiteley.com.au or call the Product Support Hotline on **1800 833 566** or New Zealand www.whiteley.co.nz Product Support Hotline **0800 257 352**

Whiteley
Medical 