

MEDICLEAN

READY TO USE HOSPITAL GRADE NEUTRAL DETERGENT

ARTG No. 125529

BENEFITS:

- Provides a streak free finish
- Neutral pH
- Non corrosive, non toxic
- Improves Infection Control by avoiding contamination through re-used containers
- Packaging 100% recyclable
-

MEDICLEAN is a ready to use neutral detergent that gives a streak free finish on all surfaces. MEDICLEAN uses Whiteley Medical's emulsion systems technology which enhances MEDICLEAN's ability to remove oils and dirt that are encountered in everyday cleaning.

MEDICLEAN is a Class 1 Medical Device and can be used to clean other Medical Device surfaces. MEDICLEAN has been formulated to replace most cleaning products and is totally compatible for cleaning before disinfection with Viraclean, hypochlorites and other hospital grade disinfectants. Do not mix with chemicals.

MEDICLEAN aids in the Infection Control process because of its' single use packaging which eliminates the need for decanting or dispensing of detergents. The elimination of Detergent dispensing systems means a reduction of risk associated with bacterial growth and contamination of Dispensing Systems and refillable bottles. Coupled with Viraclean this makes for simplification of the cleaning and disinfection process.

MEDICLEAN packaging is 100% recyclable which means it has a low environmental impact.

MEDICLEAN can be used on all surfaces including the following areas:

- Cleaning Class 1 Medical Devices
- Nursery
- CSSD
- Intensive Care Units
- Theatre Equipment
- Ward areas
- Telephones

DIRECTIONS FOR USE:

Refer to label prior to first use. Apply onto surfaces and wipe off for a streak-free finish.

Mediclean - June2011 1

CHEMICAL COMPOSITION AND PROPERTIES:

Surfactant System:.....	Blended
Colour:.....	Pale Pink
pH:.....	Neutral
Solubility:.....	100% in water
Specific Gravity:.....	1.00
Abrasives:.....	Nil
Free Caustics:.....	Nil
Phosphates:.....	Nil
Toxicity:.....	Non Toxic
Occupational Safety:.....	Safe to Use
Dangerous Goods Rating:.....	Nil
Shelf Life:.....	15 months
Safety Equipment:.....	Gloves Recommended

PACK SIZE: 12 x 500mL

HEALTH AND SAFETY STATEMENT:

There is no way of predicting an individual's response to cleaning or disinfecting products. All cleaning and disinfecting products have a potential to cause skin and other irritation in susceptible individuals. Direct physical contact with the skin should be avoided. If irritation does occur when using this product or any other cleaning/disinfecting solution then the individual should be removed from direct contact with the solution until a satisfactory safe working environment for that individual is determined.

ENVIRONMENTAL STATEMENT:

Whiteley Corporation manufactures products that comply with Australian and International environmental policy. Whiteley Corporation aims to ensure that every product is 100% biodegradable, whilst achieving optimal performance and benefit to customers. All Whiteley Corporation products are supplied in 100% recyclable containers. Whiteley Corporation and its staff are fully committed to ensuring that our Australian environment and those international environs in which our products are used are enhanced by the use of technology, and that these environments are preserved for the enjoyment of future generations.

© This document remains property of Whiteley Corporation Pty Ltd.
Alterations are not permitted without prior written authorisation from Whiteley Corporation Pty Ltd.

Mediclean - June2011 2