

LEMEX

Cleaner/Reodorant

ARTG NO. 125529



Suitable for use on



FLOORS



BENCH TOPS



NON-CRITICAL
MEDICAL DEVICES

Neutral detergent cleaners with deodorants are often recommended in hospitals for several important reasons:

1. Hygiene and Infection Control:

Neutral detergents are effective at removing dirt, oils, and microbes from surfaces without damaging the surfaces or leaving behind harmful residues. Using these cleaners helps reduce the risk of healthcare-associated infections (HAIs).

2. Non-Irritating:

Neutral pH cleaners are gentle and less likely to cause skin irritation or allergies, which is crucial in a hospital environment where patients and staff may have sensitive skin or respiratory conditions.

3. Odour Control:

Deodorising agents in cleaners help neutralise unpleasant odors, which is particularly important in medical settings.

4. Surface Protection:

Since neutral detergents are typically milder, they are less likely to damage surfaces, equipment, or finishes, which can be costly to repair or replace.

5. Regulatory Compliance:

Hospitals are often required by health authorities to use cleaning products that are both effective in disinfecting and safe for use in sensitive environments. Neutral detergent cleaners with deodorants meet these regulatory standards, ensuring compliance with health and safety guidelines.

Packsizes: Product Code:
500mL empty bottle 610536
15L 120145
2 x 5L 120146



Australian Made
Fast local supply