

# Economical, multi-enzyme detergent for complex medical instruments



**Cost effective**



**Neutral pH**



**Phosphate free**

Medizyme is a neutral pH, phosphate-free, multi-enzyme detergent designed to deliver effective reprocessing of medical devices.

Medizyme is a specialised enzymatic detergent intended for the removal of blood, protein and tissue that is found on endoscopic and surgical instruments. The proteolytic enzymes rapidly degrade proteinaceous contamination permitting easy cleaning of complex instruments.

## Medizyme is compatible for use on:

- Endoscopic equipment
- Cardiovascular instruments
- Orthopaedic instruments
- Respiratory therapy equipment
- General surgical instruments
- Dental equipment

*\*Check manufacturer's directions and warranties on complex equipment such as fiberoptical equipment before using any "wet" cleaning method.*

**Available in 500mL and 5L**



Size	Code
12 x 500mL	131034
2 x 5L	131042

ABN No. 24 000 906 678  
PO Box 1076, North Sydney NSW 2059  
Product Support Hotline: 1800 833 566  
[www.whiteley.com.au](http://www.whiteley.com.au)

NZBN: 9429041224243  
PO Box 22-519, Otahuhu, Auckland 1640  
Product Support Hotline: 0800 257 352  
[www.whiteley.co.nz](http://www.whiteley.co.nz)

Whiteley 

Protecting lives and livelihoods since 1933