

Microl

Rinsing and drying aid



Safe for use on surgical instruments



Safe for use on surfaces



70% Ethanol Solution



ARTG No: 127945

MICROL is now the benchmark standard product for device rinsing and was designed to conform with the latest guidelines. MICROL is made under strictly controlled quality control provisions in the most modern specialty Medical Devices premises in Australia.

- Filtered to 0.2µm
- Contains no dyes
- Contains 70% solution of ethanol in de-ionised water

| Size | Code |
|------------|--------|
| 12 x 500mL | 130192 |