

Product Selection Guide

Cleaning and Disinfecting Reusable Medical Devices

Choosing the right products is critical to the efficacy of every infection prevention plan.*

Cleaning



- As per manufacturer's instructions

Class 1 Medical Device

Matrix® range
Medizyme



ARTG NO: 125529

ARTG No: 125529

Disinfection

Non Critical



- Non-invasive devices
i.e. stethoscope

Low Level Instrument Grade Disinfectant

Instrumax® Pink, V-Wipes
Surfex®



ARTG No: 135544

ARTG No: 354053

ARTG No: 257360

Semi Critical



- Invasive medical device
i.e. Endoscope

High Level Instrument Grade Disinfectant

Proxy A and P**
Opal®



ARTG No: 274095

ARTG No: 206575

Sterilisation

Critical



- Invasive device that penetrates skin
i.e. Surgical instruments

Sterilants

Proxy A and P**
Aidal® Plus



ARTG No: 274095

ARTG No: 135545

Learn more or find a distributor

Australia	ABN No 24 000 906 678	New Zealand	NZBN 9429041224243
📞 1800 833 566		📞 0800 257 352	
🌐 www.whiteley.com.au		🌐 www.whiteley.co.nz	
PO Box 1076, North Sydney NSW 2059		PO Box 22-519, Otahuhu, Auckland 1640	

For further information or training please
contact your Whiteley sales representative.

*Reusable medical devices must be disinfected using
instrument grade disinfectants.

**Suitable for use in select AERs.

Whiteley



Protecting lives and livelihoods since 1933