




Water Hardness Test Kit
Product Code: 303004
\$150+GST



HACH Chloride Test Strips (40 strips)
Product Code: 303022
\$204³⁷+GST
*Suitable for testing Chloride levels between 30-600ppm.



Conductivity Meter
Product Code: 303025
\$296+GST

Endoscopy Water Quality Requirements for Compliance AS/NZS4187:2014 amd 2:2019

Water testing for Endoscopy

As a safety requirement, healthcare organisations are required to regularly test and monitor their water, in accordance with AS/NZS4187:2014 amd 2:2019.

This includes monthly water hardness, conductivity and chloride testing, and annual endotoxin testing.

Whiteley offers water testing kits and endotoxin compliance testing, to help your endoscopy and sterilising departments achieve compliance.

Whiteley also offers practical compliance support. We can train your staff and provide ongoing support, to make sure you follow industry best practice.

For testing items such as conductivity, chloride and hardness, these are to be done in accordance with the AER manufacturers recommendations.

Annual endoscopy endotoxin testing for compliance with AS/NZS4187:2014 amd 2:2019			
Final Rinse Water Endotoxin Testing*	Endotoxin 30 EU/mL ⁻¹	Annual Endotoxin Testing Fee Product Code: 303024	\$200 + GST

*Microbiological Compliance testing is also required monthly and should be conducted with your normal microbiology test provider.