




**Water Hardness Test Kit**  
Product Code: 303004  
**\$150+GST**



**HACH Chloride Test Strips (40 strips)**  
Product Code: 303022  
**\$204<sup>37</sup>+GST**  
\*Suitable for testing Chloride levels between 30-600ppm.



**Conductivity Meter**  
Product Code: 303025  
**\$296+GST**

## CSSD Water Quality Requirements for Compliance AS/NZS4187:2014 amd 2:2019

### Water testing for CSSD

As a safety requirement, healthcare organisations are required to regularly test and monitor their water, in accordance with AS/NZS4187:2014 amd 2:2019.

This includes monthly water hardness, conductivity and chloride testing, and annual endotoxin testing.

Whiteley offers water testing kits and endotoxin compliance testing, to help your endoscopy and sterilising departments achieve compliance.

Whiteley also offers practical compliance support. We can train your staff and provide ongoing support, to make sure you follow industry best practice.

AS/NZS4187:2014 amd 2:2019 WATER QUALITY USED FOR REPROCESSING OF REUSABLE MEDICAL DEVICES		
Substance	Water sample	Maximum concentration levels
<b>Water Hardness</b>	Supply Water	150mg/L
<b>Chloride**</b>	Supply Water	120mg/L
<b>Conductivity</b>	Final Rinse Water	<30µs/cm

Annual endotoxin testing for compliance with AS/NZS4187:2014 amd 2:2019			
<b>Final Rinse Water Endotoxin Testing*</b>	Endotoxin 0.25 EU/mL <sup>1</sup>	Annual Endotoxin Testing Fee Product Code: 303024	<b>\$200 + GST</b>

\*Microbiological Compliance testing is also required monthly and should be conducted with your normal microbiology test provider.

\*\*Contact your Whiteley™ Sales Representative for further information when testing below 20mg/L.