

PROXY

Automatic Endoscope Reprocessing [AER] Chemistry System*

Proudly
Australian Made



Extremely
Cost Efficient



Insist on Whiteley when ordering your Endoscopy products



ARTG No: 274095

Proxy A

Anticorrosive Additive

- For use with PROXY P.
- For use in Soluscope machines.

Size	Code
2 x 5L	160622



ARTG No: 274095

Proxy P

High Level Disinfectant
and Sterilant

- 5% peracetic acid based high level disinfectant and sterilant.
- For use with PROXY A.

Size	Code
2 x 5L	160612



ARTG No. 125529

Proxy CLN

Multi-Enzyme Detergent

- Low foaming detergent compatible with PROXY automatic endoscope reprocessing chemistry system.

Size	Code
2 x 5L	160652



ARTG No: 127945

Proxy L

70% Ethanol Solution

- 70% solution of ethanol in water.

Size	Code
4 x 2L	160632

*Suitable for use in Soluscope 3, and Sprint PAA AERs, please contact your Whiteley representative prior to first use.