

Product Selection Guide

Cleaning and Disinfecting Reusable Medical Devices

Cleaning



- As per manufacturer's instructions

Class 1 Medical Device

Clinimax®, Matrix® range,
Medizyme, Sonidet®



ARTG No: 125528 ARTG NO: 125529 ARTG No: 125529

Disinfection

Non Critical



- Non-invasive devices
i.e. stethoscope, blood pressure cuff

Low Level Instrument Grade Disinfectant

Instrumax® Pink, V-Wipes
Surfex®



ARTG No: 135544 ARTG No: 354053 ARTG No: 257360

Semi Critical



- Medical device
i.e. endoscope

High Level Instrument Grade Disinfectant

Proxy A and P
Opal®
Endox A® and Endox B®



ARTG No: 274095 ARTG No: 206575 ARTG No: 409191

Sterilisation

Critical



- Invasive device that penetrates skin
i.e. surgical instruments

Sterilants

Proxy A and P
Aidal® Plus
Endox A® and Endox B®



ARTG No: 274095 ARTG No: 135545 ARTG No: 409191

Learn more or find a distributor

Australia	ABN No 24 000 906 678	New Zealand	NZBN 9429041224243
📞 1800 833 566		📞 0800 257 352	
🌐 www.whiteley.com.au		🌐 www.whiteley.co.nz	
PO Box 1076, North Sydney NSW 2059		PO Box 22-519, Otahuhu, Auckland 1640	

For further information or training please
contact your Whiteley sales representative.

Whiteley

Protecting lives and livelihoods since 1933